

Goat KAPPA LIGHT CHAIN Polyclonal Antibody | anti-KAPPA antibody

GOAT ANTI HUMAN KAPPA LIGHT CHAIN (MOUSE ADSORBED)

CATALOG # AAA50041	HOST Goat	APPLICATIONS Western Blot, Immunofluorescence, ELISA
------------------------------	---------------------	--

PRODUCT SYNONYMS

KAPPA LIGHT CHAIN; Polyclonal Antibody; GOAT ANTI HUMAN KAPPA LIGHT CHAIN (MOUSE ADSORBED); anti-KAPPA antibody

PRODUCT SPECIFICATIONS

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Specificity	This antibody recognizes human kappa light chains. It has been cross adsorbed against pooled human myeloma proteins with lambda light chains and mouse immunoglobulins.
Form/Format	Purified Purified IgG - liquid
Concentration	IgG concentration 1.0 mg/ml (varies by lot)
Sequence Length	219
Applicable Applications for anti-KAPPA antibody	WB (Western Blot), IF (Immunofluorescence), ELISA
Immunogen	Human myeloma proteins with kappa light chains. Buffer Solution
Antiserum Preparation	Antisera to human kappa chains were raised by repeated immunisation of goats with purified antigen. Purified IgG was prepared from whole serum by affinity chromatography.
Target Species	Human
Preparation and Storage	Store at -20 degree C only. Storage in frost-free freezers is not recommended. This product should be stored undiluted. Avoid repeated freezing and thawing as this may denature the antibody. Should this product contain a precipitate we recommend microcentrifugation before use. Shelf Life: 18 months from date of despatch.

ADDITIONAL INFORMATION

Related Product Information for anti-KAPPA antibody	This antibody recognizes human kappa light chains. It has been cross adsorbed against pooled human myeloma proteins with lambda light chains and mouse immunoglobulins.
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	416539
UniProt Accession #	A2P1G9
Molecular Weight	24,086 Da
NCBI Official Full Name	kappa light chain, partial
UniProt Protein Name	Kappa light chain

UniProt Gene Name	Igkc
UniProt Entry Name	A2P1G9_MOUSE

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/polyclonal/kappa-light-chain-antibody-western-blot-immunofluorescence-elisa/50041>