

Human PSA ELISA Kit | PSA elisa kit

Human PSA ELISA Kit

CATALOG #
AAA54525

REACTIVITY
Human


PRODUCT SYNONYMS

PSA; Human PSA ELISA Kit; PSA elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	261
Samples	Serum, plasma or cell culture supernatant
Assay Type	Quantitative Sandwich
Detection Range	10 ng/ml-0.156 ng/ml
Sensitivity	Up to 0.05 ng/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for PSA elisa kit

Principle of the Assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is human PSA monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	4502173
NCBI GeneID	354
NCBI Accession #	NP_001639.1

NCBI GenBank Nucleotide #	NM_001648.2
Molecular Weight	25,116 Da
NCBI Official Full Name	prostate-specific antigen isoform 1 preproprotein
NCBI Official Synonym Full Names	kallikrein-related peptidase 3
NCBI Official Symbol	KLK3
NCBI Official Synonym Symbols	APS; PSA; hK3; KLK2A1
NCBI Protein Information	prostate-specific antigen; seminin; P-30 antigen; kallikrein-3; semenogelase; gamma-seminoprotein; prostate specific antigen
UniProt Protein Name	Prostate-specific antigen
UniProt Gene Name	KLK3
UniProt Synonym Gene Names	APS; PSA; Semin
UniProt Entry Name	KLK3_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/psa/54525>