

Mouse Myeloperoxidase (MPO) ELISA Kit | MPO elisa kit

Mouse Myeloperoxidase (MPO) ELISA Kit

CATALOG #
AAA54957

REACTIVITY
Mouse


PRODUCT SYNONYMS

Myeloperoxidase (MPO); Mouse Myeloperoxidase (MPO) ELISA Kit; MPO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	No cross-reaction with other factors.
Sequence Length	745
Samples	Mouse serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant MPO concentration
Assay Type	Sandwich
Detection Range	2000 pg/ml-31.2 pg/ml
Sensitivity	12 pg/ml.
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for MPO elisa kit

Principle of the assay: This experiment use double-sandwich elisa technique and the ELISA Kit provided is typical. The pre-coated antibody is Mouse MPO monoclonal antibody and the detecting antibody is polyclonal antibody with biotin labeled. Samples and biotin labeling antibody are added into ELISA plate wells and washed out with PBS or TBS. Then Avidin-peroxidase conjugates are added to ELISA wells in order; Use TMB substrate for coloring after reactant thoroughly washed out by PBS or TBS. TMB turns into blue in peroxidase catalytic and finally turns into yellow under the action of acid. The color depth and the testing factors in samples are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #

[188658](#)

NCBI GeneID	4353
Molecular Weight	87,248 Da
NCBI Official Full Name	myeloperoxidase
NCBI Official Synonym Full Names	myeloperoxidase
NCBI Official Symbol	MPO
NCBI Protein Information	myeloperoxidase
UniProt Protein Name	Myeloperoxidase
UniProt Gene Name	MPO
UniProt Synonym Gene Names	MPO
UniProt Entry Name	PERM_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/myeloperoxidase-mpo/54957>