

Human Complement Component 4 (C4) ELISA Kit | C4 elisa kit

Human Complement Component 4 (C4) ELISA Kit

CATALOG #
AAA54992

REACTIVITY
Human

PRODUCT SYNONYMS

Complement Component 4 (C4); Human Complement Component 4 (C4) ELISA Kit; C4 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	No cross-reaction with other factors.
Sequence Length	597
Assay Type	Sandwich
Detection Range	100 ng/mL-1.56 ng/mL
Sensitivity	Up to 0.5 ng/mL
Intra-assay Precision	<= 8%
Inter-assay Precision	<= 12%
Preparation and Storage	Store all reagents at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for C4 elisa kit

Principle of the Assay: As mentioned above, this kit utilizes the Double Antibody Sandwich ELISA technique. The pre-coated antibody is an anti-Human C4 monoclonal antibody, while the detection antibody is a biotinylated polyclonal antibody. Samples and biotinylated antibodies are added into ELISA plate wells and washed out with PBS or TBS after their respective additions to the wells. Then Avidin-peroxidase conjugates are added to the wells in after. TMB substrate is used for coloration after the enzyme conjugate has already been thoroughly washed out of the wells by PBS or TBS. TMB reacts to form a blue product from the peroxidase activity, and finally turns to yellow after addition of the stop solution (Color Reagent C). The color intensity and quantity of target analyte in the sample are positively correlated.

NCBI AND UNIPROT INFORMATION

NCBI GI #	18490166
NCBI GeneID	722
Molecular Weight	67,033 Da
NCBI Official Full Name	Complement component 4 binding protein, alpha
NCBI Official Synonym Full Names	complement component 4 binding protein, alpha
NCBI Official Symbol	C4BPA
NCBI Official Synonym Symbols	PRP; C4BP
NCBI Protein Information	C4b-binding protein alpha chain; proline-rich protein
UniProt Protein Name	C4b-binding protein alpha chain

UniProt Gene Name	C4BPA
UniProt Synonym Gene Names	C4BP; C4bp; PRP
UniProt Entry Name	C4BPA_HUMAN

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/complement-component-4-c4/54992>