

Human Glycated hemoglobin A1c ELISA Kit | GHbA1c elisa kit

Human Glycated hemoglobin A1c (GHbA1c) ELISA Kit

CATALOG #
AAA55407

REACTIVITY
Human

PRODUCT SYNONYMS

Glycated hemoglobin A1c; Human Glycated hemoglobin A1c (GHbA1c) ELISA Kit; GHbA1c elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human GHbA1c. No significant cross-reactivity or interference between Human GHbA1c and analogues was observed.
Samples	Serum, Plasma, Other biological fluids
Assay Type	Sandwich
Sample Volume	1-200 uL
Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV (%) = SD/meanX100 Intra-Assay: CV Inter-Assay: CV
Detection Wavelength	450 nm
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for GHbA1c elisa kit

Principle of the Assay: This assay employs a two-site sandwich ELISA to quantitate GHbA1c in samples. An antibody specific for GHbA1c has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyGHbA1c present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for GHbA1c is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of GHbA1c bound in the initial step. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/ghba1c-human-elisa-kits/glycated-hemoglobin-a1c/55407>