

Fish Complement 3 ELISA Kit | C3 elisa kit

Fish Complement 3 (C3) ELISA Kit

CATALOG # AAA55459	REACTIVITY Fish
------------------------------	---------------------------

PRODUCT SYNONYMS

Complement 3; Fish Complement 3 (C3) ELISA Kit; C3 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Specificity	This assay has high sensitivity and excellent specificity for detection of Fish KT. No significant cross-reactivity or interference between Fish KT and analogues was observed.
Samples	Serum, Plasma And Other Biological Fluids
Assay Type	Quantitative Competitive
Detection Range	36 pg/mL-1800 pg/mL
Sensitivity	36 pg/mL (mean of 6 independent assays).
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV(%) = SD/meanX100. Inter-Assay: CV<10%
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for C3 elisa kit	Intended Uses: For the quantitative detection of Fish 11-keto Testosterone (KT) concentration in serum, plasma and other biological fluids. Principle of the Assay: This assay employs the competitive enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with an antibody specific for KT. Standards or samples are added to the appropriate microtiter plate wells with Horseradish Peroxidase (HRP) conjugated KT. The competitive inhibition reaction is launched between with HRP labeled KT and unlabeled KT with the antibody. A substrate solution is added to the wells and the color develops in opposite to the amount of KT in the sample. The color development is stopped and the intensity of the color is measured.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/c3-fish-elisa-kits/complement-3/55459>