

Fish Epinephrine/adrenaline ELISA Kit | EPI elisa kit

Fish Epinephrine/adrenaline (EPI) ELISA Kit

CATALOG #
AAA55555

REACTIVITY
Fish

PRODUCT SYNONYMS

Epinephrine/adrenaline; Fish Epinephrine/adrenaline (EPI) ELISA Kit; EPI elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Fish
Specificity	This assay has high sensitivity and excellent specificity for detection of Fish EPI. No significant cross-reactivity or interference between Fish EPI and analogues was observed.
Assay Type	Quantitative Competitive
Samples	Serum, plasma and other biological fluids
Detection Range	12.5 pg/mL - 200 pg/mL
Sensitivity	10 pg/mL (mean of 6 independent assays)
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for EPI elisa kit

Principle of the Assay: This assay employs the competitive inhibition enzyme immunoassay technique. The microtiter plate provided in this kit has been pre-coated with an antibody specific to EPI. Standards or samples are added to the appropriate microtiter plate wells with biotin-conjugated EPI. A competitive inhibition reaction is launched between EPI (Standards or samples) and biotin-conjugated EPI with the pre-coated antibody specific for EPI. The more amount of EPI in samples, the less antibody bound by biotin-conjugated EPI. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Substrate solution is added to the wells and the color develops in opposite to the amount of EPI in the sample. The color development is stopped and the intensity of the color is measured.

Intended Uses: For the quantitative detection of Fish Epinephrine/adrenaline (EPI) concentration in serum, plasma and other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/epi-fish-elisa-kits/epinephrine-adrenaline/55555>