

Mouse Carboxyhemoglobin ELISA Kit | HbCO elisa kit

Mouse Carboxyhemoglobin (HbCO) ELISA Kit

CATALOG #
AAA55707

REACTIVITY
Mouse

PRODUCT SYNONYMS

Carboxyhemoglobin; Mouse Carboxyhemoglobin (HbCO) ELISA Kit; HbCO elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of Mouse BALP. No significant cross-reactivity or interference between Mouse BALP and analogues was observed.
Samples	Serum, Plasma And Other Biological Fluids
Assay Type	Quantitative Sandwich
Detection Range	15.6 pg/mL-1000 pg/mL
Sensitivity	7.8 pg/mL (mean of 6 independent assays).
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV (%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for HbCO elisa kit

Intended Uses: For the quantitative detection of Mouse Bone-specific alkaline phosphatase (BALP) concentration in serum, plasma and other biological fluids.

Principle of the Assay: This assay employs a two-site sandwich ELISA to quantitate BALP in Mouse serum, plasma. An antibody specific for BALP has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any BALP present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for BALP is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of BALP bound in the initial step. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/hbco-mouse-elisa-kits/carboxyhemoglobin/55707>