

Human Lyso-platelet activating factor ELISA Kit | lyso-PAF elisa kit

Human Lyso-platelet activating factor (lyso-PAF) ELISA Kit

CATALOG #
AAA55794

REACTIVITY
Human

PRODUCT SYNONYMS

Lyso-platelet activating factor; Human Lyso-platelet activating factor (lyso-PAF) ELISA Kit; lyso-PAF elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human LXA4. No significant cross-reactivity or interference between Human LXA4 and analogues was observed.
Samples	Serum, Plasma And Other Biological Fluids
Assay Type	Quantitative Sandwich
Detection Range	0.156 ng/mL-10 ng/mL
Sensitivity	0.052 ng/mL (mean of 6 independent assays).
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV (%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for lyso-PAF elisa kit

Intended Uses: For the quantitative detection of Human Lipoxin A4 (LXA4) concentration in serum, plasma and other biological fluids.

Principle of the Assay: This assay employs a two-site sandwich ELISA to quantitate LXA4 in Human serum, plasma. An antibody specific for LXA4 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any LXA4 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for LXA4 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of LXA4 bound in the initial step. The color development is stopped and the intensity of the color is measured.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/lyso-paf-human-elisa-kits/lyso-platelet-activating-factor/55794>