

Human Resolvin D1 (RvD1) ELISA Kit | RvD1 elisa kit

Human Resolvin D1 (RvD1) ELISA Kit

CATALOG #
AAA55947

REACTIVITY
Human

PRODUCT SYNONYMS

Resolvin D1 (RvD1); Human Resolvin D1 (RvD1) ELISA Kit; RvD1 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human RvD1. No significant cross-reactivity or interference between Human RvD1 and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and other biological fluids
Detection Range	62.5 pg/mL - 4000 pg/mL
Sensitivity	20.8 pg/mL (mean of 6 independent assays)
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV (%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	Store at 2-8 degree C.

ADDITIONAL INFORMATION

Related Product Information for RvD1 elisa kit

Principle of the Assay: This assay employs a two-site sandwich ELISA to quantitate RvD1 in Human serum, plasma. An antibody specific for RvD1 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any RvD1 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for RvD1 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of RvD1 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Intended Uses: For the quantitative detection of Human Resolvin D1 (RvD1) concentration in serum, plasma and other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/resolvin-d1-rvd1/55947>