

Human Interferon alpha (IFN-alpha) ELISA Kit | IFN-alpha elisa kit

Human Interferon alpha (IFN-alpha) ELISA Kit

CATALOG #
AAA56008

REACTIVITY
Human

PRODUCT SYNONYMS

Interferon alpha (IFN-alpha); Human Interferon alpha (IFN-alpha) ELISA Kit; IFN-alpha elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Specificity	This assay has high sensitivity and excellent specificity for detection of Human IFN-alpha. No significant cross-reactivity or interference between Human IFN-alpha and analogues was observed.
Assay Type	Quantitative Sandwich
Samples	Serum, plasma and other biological fluids
Detection Range	12.5 pg/mL - 800 pg/mL
Sensitivity	6.25 pg/mL (mean of 6 independent assays)
Intra-assay Precision	Intra-assay Precision (Precision within an assay) Three samples of known concentration were tested twenty times on one plate to assess intra-assay precision. Intra-Assay: CV<8%
Inter-assay Precision	Inter-assay Precision (Precision between assays) Three samples of known concentration were tested in forty separate assays to assess inter-assay precision. CV (%) = SD/meanX100. Inter-Assay: CV<12%

ADDITIONAL INFORMATION

Related Product Information for IFN-alpha elisa kit

Principle of the Assay: This assay employs a two-site sandwich ELISA to quantitate IFN-alpha in Human serum, plasma. An antibody specific for IFN-alpha has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and any IFN-alpha present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for IFN-alpha is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of IFN-alpha bound in the initial step. The color development is stopped and the intensity of the color is measured.

Intended Uses: For the quantitative detection of Human Interferon alpha (IFN-alpha) concentration in serum, plasma and other biological fluids.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/interferon-alpha-ifn-alpha/56008>