

Clostridium difficile Toxin A Monoclonal Antibody

MAb TO C. DIFFICILE TOXIN A

CATALOG # AAA58010	HOST Host: Mouse; Source: Ascites	APPLICATIONS ELISA
------------------------------	---	------------------------------

PRODUCT SYNONYMS

Clostridium difficile Toxin A; Monoclonal Antibody; MAb TO C. DIFFICILE TOXIN A; anti-Clostridium difficile Toxin A antibody

PRODUCT SPECIFICATIONS

Host	Host: Mouse; Source: Ascites
Clonality	Monoclonal
Isotype	IgG1
Specificity	Clostridium difficile Toxin A Specific for Clostridium difficile Toxin A. Does not cross react with Clostridium difficile Toxin B.
Purity/Purification	> 90% pure. Protein A Chromatography
Form/Format	Purified, Liquid
Concentration	2.1 mg/mL (OD280nm, E ^{0.1%} =1.3) (varies by lot)
Applicable Applications for anti-Clostridium difficile Toxin A antibody	ELISA
Note	Centrifuge before opening to ensure complete recovery of vial contents.
Immunogen	Purified Clostridium difficile Toxoid A.
Affinity Constant	Not Determined
Buffer	0.01 M PBS, pH 7.2 Product contains no stabilizing proteins.
Preservative	0.1% Sodium Azide
Warning	This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 - 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.
Preparation and Storage	Store at 2-8 degree C. Centrifuge before opening to ensure complete recovery of vial contents.

ADDITIONAL INFORMATION

Related Product Information for anti-Clostridium difficile Toxin A antibody	MAb to Clostridium difficile Toxin A Monoclonal Antibody to Clostridium difficile Toxin A
Product Categories/Family for anti-Clostridium difficile Toxin A antibody	Monoclonal Antibodies to Infectious Agents and Toxins

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.