

Rabbit Fluorescein Isothiocyanate (FITC) Antibody

Rabbit anti Fluorescein

CATALOG # AAA58181	HOST Rabbit	APPLICATIONS Western Blot, ELISA
------------------------------	-----------------------	--

PRODUCT SYNONYMS

Fluorescein Isothiocyanate (FITC); Rabbit anti Fluorescein; Rabbit Antibody to Fluorescein; anti-Fluorescein Isothiocyanate (FITC) antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Specificity	Recognizes Fluorescein. No cross-reactivity with immunoglobulin of any species.
Purity/Purification	Solid phase hapten-antibody capture. Affinity purified and absorbed antibody covalently coupled to a highly purified preparation of Horseradish Peroxidase (RZ>3). Care is taken to ensure adequate conjugation while preserving maximum enzyme activity. Free
Form/Format	HRP, Liquid
Concentration	1-2 mg/ml (OD280nm, E ^{0.1%} = 1.4) (varies by lot)
Applicable Applications for anti-Fluorescein Isothiocyanate (FITC) antibody	WB (Western Blot), ELISA
Immunogen	FITC-KLH conjugate
Buffer	PBS containing 10mg/ml BSA
Preservative	0.02% Thimerosal
Important Note	Centrifuge before opening to ensure complete recovery of vial contents.
Warnings	Use of sodium azide as a preservative will substantially inhibit the enzyme activity of horseradish peroxidase.
Preparation and Storage	Short-term (up to 6 months) store at 2 to 8°C. Long term, aliquot and store at -20°C. Avoid multiple freeze/thaw cycles.

ADDITIONAL INFORMATION

Related Product Information for anti-Fluorescein Isothiocyanate (FITC) antibody	Rabbit anti Fluorescein Rabbit Antibody to Fluorescein Horseradish Peroxidase Conjugated
Product Categories/Family for anti-Fluorescein Isothiocyanate (FITC) antibody	Polyclonal Antibodies to Assay Detection and Miscellaneous Reagents

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.