

Goat Kappa Light Chain, Free & Bound Polyclonal Antibody | anti-LOC100037688 antibody

Goat anti Kappa Light Chain

CATALOG # AAA58186	HOST Goat
------------------------------	---------------------

PRODUCT SYNONYMS

Kappa Light Chain; Free & Bound; Polyclonal Antibody; Goat anti Kappa Light Chain; Goat Antibody to Kappa Light Chain (Free & Bound); anti-LOC100037688 antibody

PRODUCT SPECIFICATIONS

Host	Goat
Clonality	Polyclonal
Specificity	Kappa Light Chain, Free & Bound Monospecific by immunoelectrophoresis when tested against pooled human plasma and twice concentrated pooled human serum.
Purity/Purification	Delipidation and ion exchange chromatography Product is 0.2um filtered.
Form/Format	Monospecific, Liquid
Concentration	Total protein: 9.2g/dL (Biuret (Cobas-Mira) Titer: 9mg/ml (Becker method) (varies by lot)
Immunogen	Purified human kappa light chain
Buffer	Tris-saline, pH 7.4
Preservative	0.1% Sodium azide
Important Note	Centrifuge before opening to ensure complete recovery of vial contents.
Warnings	This product contains sodium azide, which has been classified as Xn (Harmful), in European Directive 67/548/EEC in the concentration range of 0.1 - 1.0 %. When disposing of this reagent through lead or copper plumbing, flush with copious volumes of water to prevent azide build-up in drains.
Preparation and Storage	Short-term store at 2 to 8 degree C. Long term store at -20 degree C. Avoid multiple freeze/thaw cycles.

ADDITIONAL INFORMATION

Product Categories/Family for anti-LOC100037688 antibody	Polyclonal Antibodies to Human Immunoglobulins
--	--

NCBI AND UNIPROT INFORMATION

NCBI GI #	147904806
NCBI GeneID	100037688
NCBI Accession #	NP_001087252.1
NCBI GenBank Nucleotide #	NM_001093783.1
NCBI Official Full Name	kappa light chain

NCBI Official Symbol	LOC100037688
NCBI Protein Information	kappa light chain

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/polyclonal/kappa-light-chain-free-bound-antibody-elisa/58186>