

Rabbit Staphylococcus aureus Enterotoxin B (SEB) Polyclonal Antibody | anti-seb antibody

Rabbit A' Staph. Enterotoxin B

CATALOG # AAA58188	HOST Rabbit	APPLICATIONS ELISA
------------------------------	-----------------------	------------------------------

PRODUCT SYNONYMS

Staphylococcus aureus Enterotoxin B (SEB); Polyclonal Antibody; Rabbit A' Staph. Enterotoxin B; Rabbit Antibody to Staphylococcus aureus Enterotoxin B (SEB); anti-seb antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Clonality	Polyclonal
Specificity	Staphylococcus aureus Enterotoxin B (SEB) Specific to Staphylococcus aureus Enterotoxin B. Minimal cross-reactivity with staphylococcal enterotoxins A, C through E, ET, TSST and alpha hemolysin.
Purity/Purification	Toxin specific immunoaffinity column
Form/Format	Affinity Purified, Lyophilized Reconstitute with 1ml deionized water.
Concentration	10 ug IgG/ml - prior to lyophilization (varies by lot)
Applicable Applications for anti-seb antibody	ELISA
Immunogen	Staphylococcus aureus Enterotoxin B (SEB)
Buffer	Lyophilized from 0.01M phosphate, pH 7.4 containing 0.9% Sodium chloride
Preservative	None
Important Note	Centrifuge before opening to ensure complete recovery of vial contents.
Preparation and Storage	Store lyophilized product at 4 degree C. After reconstitution, aliquot and store at -80 degree C. If aliquoted for long term storage, fill volume should be equal to or greater than 50% of the nominal fill volume of the vial used. Avoid repeated freeze/thaw cycles

ADDITIONAL INFORMATION

Product Categories/Family for anti-seb antibody	Polyclonal Antibodies to Infectious Agents and Toxins
---	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	83308249
UniProt Accession #	Q2WBR3
NCBI Official Full Name	staphylococcal enterotoxin B (SEB)
UniProt Protein Name	Staphylococcal enterotoxin B (SEB)

UniProt Gene Name	seb
UniProt Entry Name	Q2WBR3_STAAU

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/polyclonal/staphylococcus-aureus-enterotoxin-b-seb-antibody-elisa/58188>