

Rabbit Campylobacter sp. Polyclonal Antibody

Campylobacter jejuni

CATALOG # AAA58529	HOST Rabbit	APPLICATIONS ELISA
------------------------------	-----------------------	------------------------------

PRODUCT SYNONYMS

Campylobacter sp.; Polyclonal Antibody; Campylobacter jejuni; Campylobacter jejuni Antiserum Products; anti-Campylobacter sp. antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Clonality	Polyclonal
Specificity	Whole cells
Form/Format	HRP
Applicable Applications for anti-Campylobacter sp. antibody	ELISA
Immunogen	strain 29428
Other Note	This antiserum has not been absorbed & may react with related microorganisms
Antiserum Characterization	(neat)

ADDITIONAL INFORMATION

Related Product Information for anti-Campylobacter sp. antibody	This product consists of purified IgG fraction of the above antiserum covalently coupled to a highly purified preparation of horseradish peroxidase (RZ >3) by methods developed by ViroStat. Care is taken to ensure adequate conjugation while preserving maximum enzyme activity. Free enzyme is absent. Immunoglobulin concentration in this product is 1-2 mg/ml. Estimated molar HRP:IgG substitution is 2-3. This preparation is supplied in a phosphate saline buffer stabilized with bovine serum albumin (10 mg/ml) and preserved with 0.002% thimerosal. This product is suitable for immunohistochemical, immunoblotting and EIA applications. Working strength must be determined by the user but a starting range of 1:20 - 1:200 is recommended for immunochemistry and blotting while 1:200 - 1:1,000 is recommended for enzyme immunoassays. Recommended short term (<6 months) storage is liquid at 2-8 degree C. For longer term storage, aliquot and freeze. Unit size is 1.0 ml.
Product Categories/Family for anti-Campylobacter sp. antibody	OMNITOPE Polyclonal Antibodies

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.