

Rabbit SS18-SSX Monoclonal Antibody | anti-SS18-SSX antibody

SS18-SSX (RBT-SS18-SSX), RMab, 0.1mL Concentrate

CATALOG # AAA59477	HOST Rabbit	REACTIVITY Paraffin, Frozen
------------------------------	-----------------------	---------------------------------------

PRODUCT SYNONYMS

SS18-SSX; Monoclonal Antibody; SS18-SSX (RBT-SS18-SSX); RMab; 0.1mL Concentrate; (RBT-SS18-SSX); ss18-ssx; anti-ss19-ssx; anti ss18-ssx; SSX1; SSX2B; SSX3; SSX4B; SSX5; SSX7; SSX8; SSX9; SSX2; SSX4; SSX6; cancer testis antigen 5 1; cancer testis antigen 5 2; Cancer testis antigen 5 3; Cancer testis antigen 5 4; cancer testis antigen family 5 member 1; cancer testis antigen family 5 member 2b; CT5 1; CT5 2; CT5 2b; CT5 3; CT5 4; HOM MEL 40; Protein SSX1; protein SSX2; Protein SSX3; Protein SSX4; protein SSX5; Protein SSX7; Protein SSX8; Protein SSX9; sarcoma synovial X chromosome related 1; SSRC; SSX; Synovial sarcoma X breakpoint 1; Synovial sarcoma X breakpoint 2; synovial sarcoma X breakpoint 2B; synovial sarcoma X breakpoint 3; synovial sarcoma X breakpoint 3 isoform a; synovial sarcoma X breakpoint 4B; synovial sarcoma X breakpoint 5; synovial sarcoma X breakpoint 7; synovial sarcoma X breakpoint 8; synovial sarcoma X breakpoint 9; tumor antigen HOM MEL 40; anti-SS18-SSX antibody


PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	Paraffin, Frozen
Clonality	Monoclonal
Isotype	IgG
Clone Number	RBT-SS18-SSX
Control	Synovial Sarcoma with the SS18-SSX fusion

IMAGES

Antibody Type	Dilution	Volume/QTY
Prediluted	Ready-To-Use	3.0 ml
Prediluted	Ready-To-Use	7.0 ml
Prediluted	Ready-To-Use	15.0 ml
Concentrated	1:50-1:200	0.1 ml
Concentrated	1:50-1:200	0.5 ml
Concentrated	1:50-1:200	1.0 ml
Control Slides		5

Dilution Info



IHC (Immunohistochemistry)
SS18-SSX Antibody Immunohistochemistry on an FFPE Smooth Muscle Tissue

ADDITIONAL INFORMATION

Product Categories/Family for anti-SS18-SSX antibody	Nuclear
--	---------

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.