

Mouse Indoleamine-2,3-Dioxygenase 2 (IDO2) ELISA Kit | IDO2 elisa kit

Mouse Indoleamine-2,3-Dioxygenase 2 (IDO2) ELISA Kit

CATALOG #
AAA65650

REACTIVITY
Mouse


PRODUCT SYNONYMS

Indoleamine-2; 3-Dioxygenase 2 (IDO2); Mouse Indoleamine-2; 3-Dioxygenase 2 (IDO2) ELISA Kit; INDOL1 Indoleamine-Pyrrole 2; 3 Dioxygenase-Like 1 Indoleamine 2; 3-dioxygenase-like protein 1 Indoleamine-pyrrole 2; 3-dioxygenase-like protein 1; IDO2 elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Mouse
Specificity	This assay has high sensitivity and excellent specificity for detection of IDO2. No significant cross-reactivity or interference between IDO2 and analogues was observed.
Sequence Length	405
Samples	Serum, plasma, tissue homogenates and other biological fluids.
Assay Type	Sandwich
Detection Range	0.781-50ng/mL
Sensitivity	0.32ng/mL
Intra-assay Precision	Intra-assay Precision (Precision within an assay): 3 samples with low, middle and high level IDO2 were tested 20 times on one plate, respectively. Intra-Assay: CV<10%
Inter-assay Precision	Inter-assay Precision (Precision between assays): 3 samples with low, middle and high level IDO2 were tested on 3 different plates, 8 replicates in each plate. CV(%) = SD/meanX100. Inter-Assay: CV<12%
Preparation and Storage	For unopened kit, all reagents should be kept according to the labels on vials. The TMB Substrate, Wash Buffer, Stop Solution should be stored at 4 degree C. All others should be stored at -20 degree C.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for IDO2 elisa kit	Intended Uses: The kit is a sandwich enzyme immunoassay for in vitro quantitative measurement of IDO2 in mouse serum, plasma, tissue homogenates and other biological fluids. Principle of the Assay: The microtiter plate provided in this kit has been pre-coated with an antibody specific to IDO2. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to IDO2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain IDO2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm +/- 10nm. The concentration of IDO2 in the samples is then determined by comparing the O.D. of the samples to the standard curve.
--	---

NCBI AND UNIPROT INFORMATION

NCBI GI #	347582671
NCBI GeneID	209176
NCBI Accession #	NP_666061.3
NCBI GenBank Nucleotide #	NM_145949.2
UniProt Accession #	Q8R0V5; A4UHF3; E9QKA9
Molecular Weight	45,255 Da
NCBI Official Full Name	indoleamine 2,3-dioxygenase 2
NCBI Official Synonym Full Names	indoleamine 2,3-dioxygenase 2
NCBI Official Symbol	Ido2
NCBI Official Synonym Symbols	Ido-2; Indol1; AI265623; C230043N17Rik
NCBI Protein Information	indoleamine 2,3-dioxygenase 2
UniProt Protein Name	Indoleamine 2,3-dioxygenase 2
UniProt Gene Name	Ido2
UniProt Synonym Gene Names	Indol1; IDO-2

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/mouse-elisa-kits/indoleamine-23-dioxygenase-2-ido2/65650>