

Human Fatty Acid Binding Protein 6, Ileal (FABP6) RTU ELISA Kit | FABP6 rtu elisa kit

Human Fatty Acid Binding Protein 6, Ileal (FABP6) ELISA Kit

CATALOG #
AAA70525

REACTIVITY
Human


PRODUCT SYNONYMS

Fatty Acid Binding Protein 6; Ileal (FABP6); Human Fatty Acid Binding Protein 6; Ileal (FABP6) ELISA Kit; GT; iBAP; I-BABP; I-BALB; I-BAP; ILBP; ILBP3; ILLBP; I-15P; Gastrotropin; Ileal Lipid-Binding Protein; Ileal Bile Acid Binding Protein; Intestinal 15 kDa protein; FABP6 rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	177
Assay Type	Sandwich
Samples	Serum, Plasma, Tissue Homogenates or Other Biological Fluids.
Detection Range	0.156-10ng/mL
Sensitivity	0.055ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for FABP6 rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to FABP6. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to FABP6. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain FABP6, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of FABP6 in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #	94721243
NCBI GeneID	2172
NCBI Accession #	NP_001035532.1
NCBI GenBank Nucleotide #	NM_001040442.1
UniProt Accession #	P51161; Q07DR7; Q8TBI3; Q9UGI7
Molecular Weight	19,874 Da
NCBI Official Full Name	gastrotropin isoform 1
NCBI Official Synonym Full Names	fatty acid binding protein 6
NCBI Official Symbol	FABP6
NCBI Official Synonym Symbols	ILBP; I-15P; I-BAP; ILBP3; ILLBP; I-BABP; I-BALB
NCBI Protein Information	gastrotropin
UniProt Protein Name	Gastrotropin
UniProt Gene Name	FABP6
UniProt Synonym Gene Names	ILBP; ILLBP; GT; ILBP; I-15P; I-BABP

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/human-rtu-elisa-kits/fatty-acid-binding-protein-6-ileal-fabp6/70525>