

Human Glial Fibrillary Acidic Protein (GFAP) RTU ELISA Kit | GFAP rtu elisa kit

Human Glial Fibrillary Acidic Protein (GFAP) ELISA Kit

CATALOG #
AAA70878

REACTIVITY
Human


PRODUCT SYNONYMS

Glial Fibrillary Acidic Protein (GFAP); Human Glial Fibrillary Acidic Protein (GFAP) ELISA Kit; Intermediate Filament Protein; GFAP rtu elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Sequence Length	431
Assay Type	Sandwich
Samples	Serum, Plasma, Tissue Homogenates, Cell Lysates, Cell Culture Supernates or Other Biological Fluids.
Detection Range	0.312-20ng/mL
Sensitivity	0.134ng/mL
RTU ELISA Kit's Advantage	The Ready-To-Use (RTU) Kits require fewer steps in setup and take less time to run. Detection Solutions are pre-diluted and total incubation time is shorter than many other ELISA Kits.
Preparation and Storage	For unopened kits: All the reagents should be kept at 4 degree C upon receipt. For opened kits: Once the kit is opened, the remaining reagents still need to be stored according to the above storage conditions. In addition, return the unused wells to the foil pouch containing the desiccant pack and reseal along entire edge of zip-seal.

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for GFAP rtu elisa kit

Test Principle: The microtiter plate provided in this kit has been pre-coated with an antibody specific to GFAP. Standards or samples are then added to the appropriate microtiter plate wells with a biotin-conjugated antibody preparation specific to GFAP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain GFAP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm $\pm$ 10nm. The concentration of GFAP in the samples is then determined by comparing the O.D. of the samples to the standard curve.

NCBI AND UNIPROT INFORMATION

NCBI GI #

196115290

NCBI GeneID	2670
NCBI Accession #	NP_001124491.1
NCBI GenBank Nucleotide #	NM_001131019.2
UniProt Accession #	P14136; Q53H98; Q5D055; Q6ZQS3; Q7Z5J6; Q7Z5J7; Q96KS4; Q96P18; Q9UFD0; B2RD44; D3DX59; E9PAX3
Molecular Weight	49,508 Da
NCBI Official Full Name	glial fibrillary acidic protein isoform 2
NCBI Official Synonym Full Names	glial fibrillary acidic protein
NCBI Official Symbol	GFAP
NCBI Official Synonym Symbols	ALXDRD
NCBI Protein Information	glial fibrillary acidic protein
UniProt Protein Name	Glial fibrillary acidic protein
UniProt Gene Name	GFAP
UniProt Synonym Gene Names	GFAP

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-rtu-elisa-kits/glial-fibrillary-acidic-protein-gfap/70878>