

Coxiella burnetii Phase II IgG Positive Control | IgG positive control

Coxiella burnetii Phase II IgG Positive Control

CATALOG #
AAA72676

PRODUCT SYNONYMS

Coxiella burnetii Phase II IgG; Coxiella burnetii Phase II IgG Positive Control; IgG positive control

PRODUCT SPECIFICATIONS

Form/Format	Human serum in protein containing phosphate buffer with Tween 20; liquid, 5 x 3 mL; negative for anti-HIV Ab, HBs-Ag (Hepatitis B-Virus surface antigen) and anti-HCV Ab; preservative: < 0.1 % sodium azide; colouring: 0.01 g/l Bromphenol blue.
Target Value	53,6 U ml
Target Range	26,9 U/mL bis/to 107 U/mL
Preparation and Storage	Store at 2-8 degree C for 1 year. DO NOT FREEZE

ADDITIONAL INFORMATION

Related Product Information for IgG positive control	These controls are positive control sera for monitoring qualitative antibody determination when using the Coxiella burnetii Phase II IgG ELISA. They are to be regarded as additional controls to the reagents supplied with the test kit. They are used to determine validity of the test runs, the precision and reliability of the method, and are particularly recommended as an aid for internal quality management in accredited laboratories.
--	--

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com
Online datasheet: <https://www.aaabiotech.com/positive-control/coxiella-burnetii-phase-ii-igg/72676>