

Rabbit anti-HEV US-u3-2e Hepatitis E Virus Polyclonal Antibody

Hepatitis E Virus antibody

CATALOG # AAA75006	HOST Rabbit	REACTIVITY HEV US-u3-2e	APPLICATIONS Western Blot
------------------------------	-----------------------	-----------------------------------	-------------------------------------

PRODUCT SYNONYMS

Hepatitis E Virus; Polyclonal Antibody; Hepatitis E Virus antibody; Polyclonal Hepatitis E Virus; Anti-Hepatitis E Virus; Hepatitis E; Hep E; HEV; anti-Hepatitis E Virus antibody

PRODUCT SPECIFICATIONS

Host	Rabbit
Reactivity	HEV US-u3-2e
Clonality	Polyclonal
Isotype	IgG
Specificity	HEV US-u3-2e
Purity/Purification	Hepatitis E Virus antibody was purified by affinity chromatography.
Form/Format	Supplied in PBS with 0.02% NaN ₃ .
Concentration	1.0 mg/ml (varies by lot)
Applicable Applications for anti-Hepatitis E Virus antibody	WB (Western Blot)
Biological Significance	Hepatitis E virus (HEV) is a positive-sense single-stranded RNA icosahedral virus with a 7.5 kilobase genome. HEV has a fecal-oral transmission route. It is one of five known hepatitis viruses: A, B, C, D, and E.
Biohazard Information	This product contains sodium azide as preservative. Although the amount of sodium azide is very small appropriate care must be taken when handling this product.
Immunogen	Hepatitis E Virus antibody was raised in rabbit using residues 627-637 [CPECRTLGLQG] of the 71 kDa HEV capsid protein as the immunogen.
Preparation and Storage	Store at 4 degree C for short term storage. Aliquot and store at -20 degree C for long term storage. Avoid repeated freeze/thaw cycles.

ADDITIONAL INFORMATION

Related Product Information for anti-Hepatitis E Virus antibody	Affinity purified Rabbit polyclonal Hepatitis E Virus antibody
Product Categories/Family for anti-Hepatitis E Virus antibody	Infectious Disease ; Purified Polyclonal Antibodies

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.