

Human Cardiolipin Total Ab ELISA Kit

Cardiolipin Total Ab ELISA Kit

CATALOG #
AAA77719

REACTIVITY
Human

PRODUCT SYNONYMS

Cardiolipin Total Ab; Cardiolipin Total Ab ELISA Kit; Cardiolipin Total Ab ELISA; Cardiolipin Total Ab elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Human
Samples	Human serum or plasma
Intended Uses	The anti-Cardiolipin antibody (aCL) ELISA Kit is intended for the detection of IgG, IgM and IgA antibody to Cardiolipin in human serum or plasma. For research use only.
Preparation and Storage	Product should be stored at 2-8 degree C. Product is stable for 24 months from the date of manufacturing.

ADDITIONAL INFORMATION

Product Categories/Family for Cardiolipin Total Ab elisa kit	Autoimmune Disorders
References	1. Barka N, Reagan K, Agopian M, Peter JB. Frequency of anticardiolipin and antiphosphatidylserine antibodies in patients with suspected antiphospholipid syndrome. Clin Chem 1995;41:5. 2. Khamashta MA, Hughes GRV. Detection and importance of anticardiolipin antibodies. J Clin Pathol 1993;46:104-7. 3. Reyes H, Dearing L, Bick RL, Shoenfeld Y, Peter JB. Laboratory diagnosis of antiphospholipid syndromes. Lab Clin North Am 1995;15:85-108. 4. Harris EN, Pierangeli S. Anticardiolipin antibodies: specificity and function. Lupus 1994;3:217-22. 5. Barka N. Phospholipid autoantibodies - phosphatidylserine. In: Peter JB, Shoenfeld Y, editors. Autoantibodies. Amsterdam: Elsevier 1996:630-4. 6. Triplett DA. Assays for detection of antiphospholipid antibodies. Lupus 1994;3:281-7. 7. Khamashta MA, Cuadrado MJ, Mujic F, et al. The management of thrombosis in the antiphospholipid antibody syndrome. N Engl J Med 1995;332:993-7. 8. Olee T, Pierangeli S, Handley H, et al. A monoclonal IgG anticardiolipin antibody from a patient with the antiphospholipid syndrome is thrombogenic in mice. Proc Natl Acad Sci 1996;93:8606-11. 9. Hughes GRV. The antiphospholipid syndrome: ten years on. Lancet 1993;342:341-4.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | contact@aaabiotech.com

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/cardiolipin-total-ab/77719>