

## Human Measles IgG ELISA Kit

Measles IgG ELISA Kit

CATALOG #  
**AAA77740**

REACTIVITY  
**Human**

### PRODUCT SYNONYMS

Measles IgG; Measles IgG ELISA Kit; Measles IgG ELISA; Measles IgG elisa kit

### PRODUCT SPECIFICATIONS

|                         |                                                                                                           |
|-------------------------|-----------------------------------------------------------------------------------------------------------|
| Reactivity              | Human                                                                                                     |
| Samples                 | Leukocyte                                                                                                 |
| Preparation and Storage | Product should be stored at 2-8 degree C. Product is stable for 24 months from the date of manufacturing. |

### ADDITIONAL INFORMATION

|                                                       |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Related Product Information for Measles IgG elisa kit | <p>Intended Uses: Leukemia MLL-AF4 Fusion Gene real time RT-PCR Kit is used for the detection of Leukemia MLL-AF4 Fusion Gene in leukocyte by using real time PCR systems.</p> <p>Principle of the Assay: The principle of the real-time detection is based on the fluorogenic 5'nuclease assay. During the PCR reaction, the DNA polymerase cleaves the probe at the 5' end and separates the reporter dye from the quencher dye only when the probe hybridizes to the target DNA. This cleavage results in the fluorescent signal generated by the cleaved reporter dye, which is monitored real-time by the PCR detection system. The PCR cycle at which an increase in the fluorescence signal is detected initially is proportional to the amount of the specific PCR product. Monitoring the fluorescence intensities in real-time allows the detection of the accumulating product without having to re-open the reaction tube after the amplification.</p> |
| Product Categories/Family for Measles IgG elisa kit   | <a href="#">Infectious Disease</a>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

**For Research Use Only.** Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.

AAA Biotech | Tel: +1 (619) 604-9440 | Fax: +1 (619) 604-9441 | [contact@aaabiotech.com](mailto:contact@aaabiotech.com)

Online datasheet: <https://www.aaabiotech.com/human-elisa-kits/measles-igg/77740>