

GFP/HEK293 cell line

GFP/HEK293 Stable Cell Line

CATALOG # AAA78569	APPLICATIONS Flow Cytometry, Functional Assay, Functional Assay
------------------------------	---


PRODUCT SYNONYMS

GFP/HEK293; GFP/HEK293 Stable Cell Line; GFP/HEK293 cell line

PRODUCT SPECIFICATIONS


Form/Format	Each vial contains 2 ~ 3 x 10 ⁶ cells in 1 ml of 90% FBS + 10% DMSO
Applicable Applications for GFP/HEK293 cell line	FCM/FACS (Flow Cytometry), FA (Functional Assay)
Preparation and Storage	Immediately upon receipt, store in liquid nitrogen. Ship: Dry ice

IMAGES



FCM/FACS (Flow Cytometry)

Fig-2: Detection of GFP in the GFP/HEK293 stable cell line through flow cytometry . Parental HEK293 cells (Green); GFP/HEK293 cells (Red).



IF (Immunofluorescence)

Fig-1: Analysis of the GFP/HEK293 stable cell line through fluorescence microscopy. Bright-field image (Left); Fluorescence image (Right).

ADDITIONAL INFORMATION

Related Product Information for GFP/HEK293 cell line	GFP/HEK293 Stable Cell Line is a stably transfected HEK293 cell line which expresses enhanced green fluorescent protein (eGFP). Sequence data: A eGFPMVSKGEELFTGVVPILVELDGDVNGHKFSVSGEGEGDATYGLTKFICTTGKLPVPWPTLVTTLTYGVCFSRYPDHMKQHDFKSA MPEGYVQ
Product Categories/Family for GFP/HEK293 cell line	Stable Cell Lines

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.