

Porcine D-Xylose ELISA Kit

Porcine D-Xylose ELISA Kit

CATALOG #
AAA97389

REACTIVITY
Porcine


PRODUCT SYNONYMS

D-Xylose; Porcine D-Xylose ELISA Kit; D-Xylose elisa kit

PRODUCT SPECIFICATIONS

Reactivity	Porcine
Samples	Serum, plasma, cell culture supernatants, body fluid and tissue homogenate
Assay Type	Quantitative Competitive
Sensitivity	1.0 ng/mL.
Preparation and Storage	Store all reagents at 2-8 degree C

IMAGES



Standard Curve (Sample)

ADDITIONAL INFORMATION

Related Product Information for D-Xylose elisa kit

Intended Uses: This DXY ELISA kit is a 1.5 hour solid-phase ELISA designed for the quantitative determination of Porcine DXY. This ELISA kit for research use only!

Principle of the Assay: DXY ELISA kit applies the competitive enzyme immunoassay technique utilizing a monoclonal anti-DXY antibody and an DXY-HRP conjugate. The assay sample and buffer are incubated together with DXY-HRP conjugate in pre-coated plate for one hour. After the incubation period, the wells are decanted and washed five times. The wells are then incubated with a substrate for HRP enzyme. The product of the enzyme-substrate reaction forms a blue colored complex. Finally, a stop solution is added to stop the reaction, which will then turn the solution yellow. The intensity of color is measured spectrophotometrically at 450nm in a microplate reader. The intensity of the color is inversely proportional to the DXY concentration since DXY from samples and DXY-HRP conjugate compete for the anti-DXY antibody binding site. Since the number of sites is limited, as more sites are occupied by DXY from the sample, fewer sites are left to bind DXY-HRP conjugate. A standard curve is plotted relating the intensity of the color (O.D.) to the concentration of standards. The DXY concentration in each sample is interpolated from this standard curve.

For Research Use Only. Not for use in diagnostic or therapeutic procedures. Small volumes of product may become entrapped in the seal of the vial during shipment and storage; if necessary, briefly centrifuge the vial on a tabletop centrifuge to dislodge any liquid in the cap.